

IMSERC User Manual for NMR

CONTENTS

Introduction.....	3
Safety.....	3
Data Management.....	3
Software	4
Default Instrument Status	4
Sample Preparation	5
FIVE steps to use IMSERC NMR	6
Common Commands/Parameters in TopSpin	6
Setting up experiments & Processing.....	6
Important parameters for acquisition.....	7
Data Acquisition:	8
Use IconNMR w/ SampleXpress	8
Interactively setup experiment with TopSpin	10
Data Analysis	14
Publication.....	14
Experimental Section.....	14
Acknowledgement.....	14
Troubleshooting	15
APPENDIX A: List of common experiments on ICONNMR.....	16
APPENDIX B: Variable Temperature Control for NMR	17
Introduction.....	17
Operation Procedure.....	17
Important notes (instrument damage will occur if not followed):.....	19
APPENDIX C: 1D Selective NOESY with Topspin	21
Reference Spectrum	21
Set up the 1D selective NOESY.....	21
Display 1D NOESY together with 1D proton spectrum	24
APPENDIX D: Proton DOSY Experiment.....	25

Introduction.....	25
Experiment	26
Processing.....	27
Additional notes	28
APPENDIX E: NMR T1 relaxation Experiment.....	29
Introduction.....	29
Setting up proton T1 Experiment.....	29
Processing.....	32
APPENDIX F: Kinetics / reaction monitoring.....	35
Introduction.....	35
Separate 1D Spectra w/ multi_zgvd	35
Pseudo 2D Mode Procedure.....	36
APPENDIX G: Evans Method.....	38
Introduction.....	38
Preparation.....	38
Experiment Setup (quick procedures).....	39
Data Analysis	39
References.....	40
Revisions.....	41

IMSERC User Manual for NMR (V2.0)

INTRODUCTION

This manual is intended to give you instructions on how to setup the routine NMR experiments with Topspin and IconNMR interfaces. Use of this instrument is allowed only by qualified users after receiving training by a staff member. Do not run this instrument without approval from IMSERC staff. Failure to do so may cause damage to the instrument, produce invalid data, and result in additional fees and/or removal of all IMSERC privileges. This short set of instructions is meant to serve as a guide for 'routine' data collection on the instrument. For custom experiments, contact a staff member. Please read this standard operating procedure and acquaint yourself with the instrument. If during the course of using the system, something happens that you do not understand, please **stop** and **get help**. In any event, be completely prepared to justify your actions. The cost of even minor repairs is considerable.

SAFETY

All users of IMSERC must review the general safety policies at <http://imserc.northwestern.edu/about-policies.html>. To become an independent user of this instrument, you must have the following safety training and certificates that are offered at <https://learn.northwestern.edu>:

- Laboratory Safety
- Strong magnetic field safety
- Personal Protective Equipment

You need the above certificates in order to be able to reserve time for this instrument on NUCore. Upon completion of the certificate, it will take an overnight to filter through the different systems and get into the files that NUCore uses. Additionally, familiarize yourself with the location of standard safety stations like eyewash and shower stations found in the west side of the NMR room. Protective eyewear is required in this room, and gloves should be removed when using the computer.

DATA MANAGEMENT

The following template is used to save your fids: /home/walkon/data/YourGroupFolder/YourNetid/SampleName. Your personal data folder is created during training. Please save data under your personal folder, which must be located under your group folder (supervisor's last name). See a staff member if you need help to create your personal folder on the instrument.

IMSERC User Manual for NMR (V2.0)

Data on the instrument is copied to 'imsercdata.northwestern.edu' in real time. You can access your NMR data through imsercdata.northwestern.edu/YourGroupFolder/nmr/NMRInstrumentName. Please follow instructions at <http://imserc.northwestern.edu/about-general-faq.html#data> for details about data access.

SOFTWARE

Offline data processing and analysis can be performed with MNOVA and Topspin. Northwestern has campus wide license for MestraNova. Please refer to Mnova installation instructions at <http://imserc.northwestern.edu/about-general-faq.html#software> to install the program on your computer. Make sure to connect to Northwestern VPN when you work off campus. You are also encouraged to download the Topspin software from Bruker at <https://www.bruker.com/service/support-upgrades/software-downloads/nmr.html> and claim your free academic license.

DEFAULT INSTRUMENT STATUS

The default interface on Ag500 and X500 is **IconNMR**. Please do not change! You do not need go through NUCore to use these two fully automated NMR instruments.

The default interface on A600, Au400, and HFCN600 is **Topspin**. Computer screen is by default deactivated. You must start your reservation through NUCore in order to turn on the computer screen. You have to login to NUCore to start your reservation before you can use. You can optionally start IconNMR if you want. Please return the instrument back to Topspin after you finish using IconNMR. Please end your reservation on NUCore after you finish.

The Hg400 (solids) is always on **Topspin**. Similar to A600/Au400/HFCN600, you have to go through NUCore to start/end your reservations. It has two probes, 4mm HX and 1.6 mm HFX. Notify the staff member in advance, if you need use a different probe than the one is currently installed.

If there is an error or problem with the instrument, please report the issue by following at least one of the steps below:

1. If you have already started your reservation using NUCore, please logoff by selecting the error reporting option and a brief description about the issue

IMSERC User Manual for NMR (V2.0)

2. If you have not started your reservation using NUCore, please report problems with the instrument at <http://imserc.northwestern.edu/contact-issue.html> add place the 'Stop' sign near the instrument computer. 'Stop' signs are located next to instrument and online at the link above
3. Email or talk to a staff member

SAMPLE PREPARATION

Please use clean, non-scratched 5mm NMR tubes to prepare your samples. Some reputable NMR tubes vendors are listed here: <http://imserc.northwestern.edu/nmr-links.html>

The sample volume should be between 0.5 to 0.7 ml, preferably 0.55 ml for best shimming results. The solution should be clear of precipitation or suspended particles by going through filtration or centrifugation.

IMSERC User Manual for NMR (V2.0)

FIVE STEPS TO USE IMSERC NMR

1. Login to NUCORE with your netid and logon to the instrument:
 - a) turn on the computer monitor at instrument,
 - b) your usage count starts
2. Login to instrument with your operator id (usually same as netid)
3. Load your sample and run your experiment
4. Logout from the instrument
5. Login to NUCORE and logout your instrument session:
 - a) turn off the computer monitor at instrument,
 - b) your usage count stops

COMMON COMMANDS/PARAMETERS IN TOPSPIN

SETTING UP EXPERIMENTS & PROCESSING

sx 10 to put #10 (could be any position) sample on the autosampler into magnet.

sx ej to put sample inside magnet back to the autosampler

rga automatically set receiver gain

zg start acquisition

tr transfer data (while acquisition is in progress)

multizg start multiple acquisitions starting from current dataset

go submit experiment to acquisition

stop abort an acquisition, losing all the FID data recorded so far

halt halt the running acquisition, saving the recorded FID data to hard disk

efp weighted Fourier Transformation for 1D dataset

apk do automatic phase correction

abs automatically optimize baseline

IMSERC User Manual for NMR (V2.0)

xfb weighted Fourier Transformation for 2D dataset

IMPORTANT PARAMETERS FOR ACQUISITION

- P1** F1 channel 90° pulse width, micro seconds
- P2** F1 channel 180° pulse width
- RG** Receiver gain
- D1** relaxation delay, 1 to 5 times T1
- 2TD** Time domain data points for F2 (direct dimension)
- 2SW** spectral width in ppm for F2 (direct dimension)
- 1TD** Time domain data points for F1 (indirect dimension)
- 1SW** spectral width in ppm for F1 (indirect dimension)
- AQ** Acquisition time in seconds
- NS** Number of scans
- DS** Number of dummy scans
- NUC1 – NUC8** Nucleus observed (1H, 13C, 31P, 19F, etc.)
- O1 – O8** Frequency offset for channel 1 – 8 in Hz
- O1P – O8P** Freq. offset for channels 1 – 8 in ppm
- SFO1 – SFO8** Freq. for channels 1 – 8 in MHz



IMSERC User Manual for NMR (V2.0)

DATA ACQUISITION:


USE ICONNMR W/ SAMPLEXPRESS

1. Login with your operator ID
2. Load your sample to SampleExpress
3. Click the Holder # where you sample is loaded.
4. Click Add
5. Fill in following fields: Name (5a), No. (5b), Solvent (5c), and Experiment (5d)
6. Change parameters if needed
7. Click Submit
8. Logout ICON by clicking on “change user”

Fig 4: Steps to setup NMR experiment with IconNMR



IMSERC User Manual for NMR (V2.0)



Notes for middle panel of fig 4 (step #2, labeled with 1-4 in black):


1. Message display
2. Displays the cassette with its 60 positions.
 - Gray positions are empty positions.
 - White positions are loaded with samples.
3. Information about Operator Access Position and Sample inside magnet.
4. Touch screen buttons to **Add** and **Remove** samples.

IMSERC User Manual for NMR (V2.0)

INTERACTIVELY SETUP EXPERIMENT WITH TOPSPIN

1. Create new dataset and setup initial parameters
2. Insert/Eject Sample
3. Lock on your selected Solvent
4. Tune/Match the Probe to the nuclei of your experiment
5. Shim
6. Load probe related parameters
7. Auto set receiver gain
8. Acquire fid
9. Process data

Fig 1: Flow chart for interactively setting up NMR experiment with TOPSPIN



Parameters setting for 1D and 2D experiments

The arrowed fields at Fig 2 (1-5) need to be filled and changed at Fig 3 (1-8) based on your needs.

Input to the “new dataset” window:

1. Sample name
2. Experiment number
3. Experiment to run
4. Solvent used
5. Directory for the dataset (your folder under your group name)
6. Change the parameters for fields 1-8 at Fig. 3 in next page based on your needs

IMSERC User Manual for NMR (V2.0)

Fig 2: Input to the “new dataset” window

The screenshot shows a window titled "Create New Dataset - new" with a close button (X) in the top right corner. The window contains the following text and controls:

Prepare for a new experiment by creating a new data set and initializing its NMR parameters according to the selected experiment type. For multi-receiver experiments several datasets are created. Please define the number of receivers in the Options.

NAME: (Callout 1)

EXPNO: (Callout 2)

PROCNO:

Use current parameters

Experiment: (Callout 3)

Options (expanded):

- Set solvent: (Callout 4)
- Execute 'getprosol'
- Keep parameters:

DIR: (Callout 5)

Show new dataset in new window

Number of additional datasets: (1,2, ...16)

TITLE:

Buttons at the bottom:

IMSERC User Manual for NMR (V2.0)

Fig 3: Important parameters might need to be changed

The screenshot shows the 'General Channel f1' parameters in the Bruker software. The parameters are organized into two sections: 'General' and 'Channel f1'. Callouts 1 through 5 highlight specific parameters that might need to be changed:

- 1:** RG (Receiver gain) is set to 203.
- 2:** SWH [Hz, ppm] is set to 9615.38 and 16.0212.
- 3:** D1 [sec] is set to 2.00000000.
- 4:** NS is set to 4.
- 5:** SFO1 [MHz] is set to 600.1678208.

Parameter	Value	Description
PULPROG	zgpr	Pulse program for acquisition
TD	32768	Time domain size
SWH [Hz, ppm]	9615.38 / 16.0212	Sweep width
AQ [sec]	1.7039360	Acquisition time
RG	203	Receiver gain
DW [µsec]	52.000	Dwell time
DE [µsec]	6.50	Pre-scan-delay
D1 [sec]	2.00000000	Relaxation delay; 1-5 * T1
d12 [sec]	0.00002000	Delay for power switching [20 usec]
DS	0	Number of dummy scans
NS	4	1 * n, total number of scans: NS * TDO
TDO	1	Number of averages in 1D
SFO1 [MHz]	600.1678208	Frequency of ch. 1
O1 [Hz, ppm]	2820.78 / 4.700	Frequency of ch. 1
NUC1	1H	Nucleus for channel 1
P1 [µsec]	12.700	F1 channel - 90 degree high power pulse
PLW1 [W, dB]	31.623 / -15.00	F1 channel - power level for pulse (default)
PLW9 [W, dB]	0.00020402 / 36.90	F1 channel - power level for presaturation

The screenshot shows the 'Experiment' and 'Width' parameters in the Bruker software. The parameters are organized into two sections: 'Experiment' and 'Width'. Callouts 6 and 7 highlight specific parameters that might need to be changed:

- 6:** NS is set to 8.
- 7:** SW [ppm] is set to 20.0264 and 165.0000.

Parameter	Value	Description
PULPROG	hsqcetgp	Current pulse program
AQ_mod	DQD	Acquisition mode
FnTYPE	traditional(planes)	nD acquisition mode for 3D etc.
FnMODE	Echo-Antiecho	Acquisition mode for 2D, 3D etc.
TD	1024 / 256	Size of fid
DS	16	Number of dummy scans
NS	8	Number of scans
TDO	1	Loop count for 'td0'
TDav	0	Average loop counter for nD experiments
SW [ppm]	20.0264 / 165.0000	Spectral width
SWH [Hz]	12019.230 / 24902.283	Spectral width
IN_F [µsec]		Increment for delay
AQ [sec]	0.0425984 / 0.0051401	Acquisition time
FIDRES [Hz]	23.475060 / 194.549088	Fid resolution
FW [Hz]	4032000.000	Filter width
RG	203	Receiver gain
DW [µsec]	41.600	Dwell time

IMSERC User Manual for NMR (V2.0)

RR 3 1 /home/walkon/data/zhang

Spectrum ProcPars **AcquPars** Title PulseProg Peaks Integrals Sample Structure Plot Fid Acqu

Probe: PA BBO 600S3 BB-H-D-05 Z BTO

Experiment Width Receiver Nucleus Durations Power Program Probe Lists

Nucleus 1

NUC1	1H	Edit...	13C	Observe nucleus
O1 [Hz]	2820.78		11318.37	Transmitter frequency offset
O1P [ppm]	4.700		75.000	Transmitter frequency offset
SFO1 [MHz]	600.1678208		150.9229277	Transmitter frequency
BF1 [MHz]	600.1650000		150.9116093	Basic transmitter frequency

8

IMSERC User Manual for NMR (V2.0)

DATA ANALYSIS

MNOVA has an extremely comprehensive manual accessed through “help”. A simplified MNOVA manual can be found at http://imserc.northwestern.edu/downloads/nmr-mnova_chemists8_simplified.pdf. The most useful Topspin manuals for processing is at <http://imserc.northwestern.edu/downloads/nmr-processing-reference.pdf>

PUBLICATION

EXPERIMENTAL SECTION

Write a paragraph similar to the ‘experimental section’ found in a publication. Details about every instrument component and methods must be provided in this paragraph. The NMR spectrometer specifications can be found in the table at “**ACKNOWLEDGEMENT**” section below. You need only the content in the parenthesis. For example, you should write “Bruker NOE 600 MHz spectrometer equipped with QCI-F cryoprobe” if you use HFCN600, DO NOT write “HFCN600” as the name of NMR spectrometer in your publication.

The important information about your NMR experiment includes but not limited to, the experiment name, spectral width, acquisition time/time domain points, window function for apodization etc.

ACKNOWLEDGEMENT

Acknowledgement info is listed under <http://imserc.northwestern.edu/about-acknowledgements.html>.

If your work used the HFCN600 please use this acknowledgement:

This work made use of the IMSERC at Northwestern University, which has received support from the NIH (1S10OD012016-01 / 1S10RR019071-01A1), Soft and Hybrid Nanotechnology Experimental (SHyNE) Resource (NSF ECCS-1542205), the State of Illinois, and the International Institute for Nanotechnology (IIN).

For other NMR instruments, please use this acknowledgement:

For This work made use of the IMSERC at Northwestern University, which has received support from the Soft and Hybrid Nanotechnology Experimental (SHyNE) Resource (NSF ECCS-1542205), the State of Illinois, and the International Institute for Nanotechnology (IIN).

IMSERC User Manual for NMR (V2.0)

Instrument	Funding Source
HFCN600 (Bruker NOE 600 MHz w/ QCI-F cryoprobe)	NIH 1S10OD012016-01 / 1S10RR019071-01A1
A600 (Bruker AVANCE III 600 MHz w/ BBFO probe)	NU / Int. Institute of Nanotechnology
Ag500 (Bruker AVANCE III 500 MHz w/ DCH cryoprobe)	NU
X500 (Bruker AVANCE HD 500 MHz w/ Prodigy probe)	Int. Institute of Nanotechnology
Hg400 (Bruker AVANCE HD 400 MHz w/ MAS solids probe)	Int. Institute of Nanotechnology
Au400 (Bruker AVANCE HD Nanobay 400 MHz w/ BBFO probe)	NSF CHE-1048773

TROUBLESHOOTING

1. If you need go to Topspin interface by exiting IconNMR on A600, Au400, and HFCN600, please login with username “walkon” (with password “GO*****”), stop automation first, then quit IconNMR.
2. If you have to reboot the computer for NMR instrument, login with account “walkon” (with password “GO*****”). Click the “Topspin” icon on desktop to start acquisition software. Leave the acquisition software open when you are done with the measurement.
3. At this time, if you experience difficulty with SampleXpress autosampler, please contact a Staff.

IMSERC User Manual for NMR (V2.0)

APPENDIX A: LIST OF COMMON EXPERIMENTS ON ICONNMR

Experiment entries	Description	Recommended setting for ~10mg material, please modify "NS" accordingly based on your sample concentration	
PROTON_icon	routine 1D proton spectrum	<p>Common parameters for all experiments:</p> <p>D1: relaxation delay, 1 to 5 times of T1, 2-5 seconds SW: spectral width in ppm for F2 (direct dimension) AQ: Acquisition time in seconds NS: number of scans O1P: offset freq for channels 1 in ppm, usually 1H O2P: offset for channel 2, usually 13C 1TD: Time domain data points for F1 (indirect dimension), aka number of increments 1SW: spectral width in ppm for F1</p>	
C13CPD_icon	1D ¹ H-decoupled ¹³ C spectrum		
C13DEPTQ135_icon	DEPTQ 135 experiment detect all Carbons - CH3/CH positive CH2/C negative		
C13IG_icon	13C with inverse gated 1H decoupling no NOE for quantitative NMR		
C13DEPT90_icon	DEPT 90 experiment only CH		
C13DEPT135_icon	DEPT 135 experiment CH3/CH positive CH2 negative - ¹³ C 1-bond correlations, all peaks positive (dept-45 analog) DEPT-135 experiment		
gCOSY_icon	Gradient selected COSY		
COSYDQF_icon	COSY with double quantum filter		
HSQC_EDIT_icon	1H-13C multiplicity edited HSQC with gradient selection		
HSQC_icon	1H-13C 1-bond correlations, all peaks positive, HSQC with gradient selection		
HSQC_EDIT_NUS_icon	1H-13C multiplicity edited HSQC with gradient selection Non Uniform Sampling w/ 25% sampling density		
HMBC_icon	1H-13C HMBC with gradient selection using 3-fold low pass filter for better 1J suppression		CNST13 = 3-12 Hz (default 8Hz for J _{2/3} CH)
HMBC_NUS_icon	1H-13C HMBC with gradient selection using 3-fold low pass filter for better 1J suppression Non Uniform Sampling w/ 50% sampling density		
TOCSY_icon	Phase sensitive 2D TOCSY experiment using MLEV-17 mixing	d9 = 30 to 120 ms	
NOESY_icon	Phase sensitive NOESY 1H-1H correlations based on proximity also for exchange	d8 = 0.1 to 1 second	
ROESY_icon	1H-1H correlations based on proximity for intermediate MW around 1600 Da	p15 = 0.1 to 0.5 second	
WATER_SUPP_icon	Solvent suppression with noesygppld sequence		

IMSERC User Manual for NMR (V2.0)

APPENDIX B: VARIABLE TEMPERATURE CONTROL FOR NMR

INTRODUCTION

This manual is to instruct NMR IMSERC users to do the Variable Temperature (VT) experiments on A600, Au400, and Hg400 (w/ RT probes) NMR spectrometers (-150 to +150 °C). In principle, one can follow similar procedure to use the HFCN600 (equipped with QCI-F cryoprobe, with much smaller VT range from -30 to 70 °C).

Depending on your target temperature, please consult with **fig 1** and **table 1, 2** to pick the right spinner, set up the gas flow rate, and the cooling power for the BCU II device.

When operating at elevated temperatures, internal probe heating is sufficient and precisely regulated with the new BSVT on Au400 and Hg400, which also controls and regulates the selected VT flow rates and other auxiliary flows. Chillers are not required as long as the operating temperature is sufficiently above the room temperature. For A600, you need manually adjustments of Shim/Flush gas. If you need temperatures below room temperature, chiller (BCU II or FTS) is required.

Always remember to redo the lock, tune, and shimming if the temperature changes more than 10 °C. The solvent boiling and/or frozen point must be considered. For high temperature, the highest temperature allowed should be at least 10 °C below the boiling point while for a low temperature experiment; the lowest temperature should be at least 5 °C above frozen point.

The real temperature and detected temperature may be slight different. Please refer to the temperature calibration curve for correction or calibrate temperature yourself. If you need a temperature below -60 °C, please ask NMR staff to use liquid nitrogen as a cooling source.

OPERATION PROCEDURE

1. Select the right spinner based on **fig 1** and **table 1**, position and load your sample.
2. Open temperature control panel by typing command **edte**.
3. To change gas flow rate, click "**set**" under "Target Gas Flow" on BCU II (**fig 2**), put appropriate gas flow rate in the popup box. For FTS chiller, simple click on "-" or "+" to change the flow rate (**fig 3**).

IMSERC User Manual for NMR (V2.0)

4. To change the chiller power, for FTS chiller, set the desired temperature using the ▲▼ arrow. For BCU II, click “set” under “Target Power” and select cooling power based on **table 2**.
5. To change temperature, click “set” under “Target Temperature” on BCU II (**fig 2**) to set your target. For FTS, click “change” button of “Target temp” row. Change the temperature by 10 °C increment and give around 5 minutes before doing next increment. Changing temperature too fast may damage the probe.
6. For FTS chiller, you can explicitly set the temperature. The maximum you can set is 60 °C. Please keep FTS at least 10 °C below your target temperature for above RT experiment. For extremely low experiment, you have to keep the FTS temperature 30 to 40 °C below your target.
7. Once the target temperature is regulated, please wait at least 5-10 minutes for the sample to reach equilibrium before you run your experiment. Remember to redo locking, tuning, and shimming for each different temperature.
8. Collect data.
9. After finishing, please restore the temperature setting back to default at 25 °C. Change flow rate and cooling power back to default.

Fig 1: Three types of spinners: Standard **POM** Spinner (blue, left), 0°C to +80°C; **Kel-F** Spinner (off white, middle) for elevated Temperatures, +80°C to +120°C; **Ceramics** Spinner (white, right) for high and low Temperatures, +120°C to +180°C and 0°C to -150°C (A600, Au400, and Hg400)



IMSERC User Manual for NMR (V2.0)

Table 1: The VT gas, shim gas, flush gas settings and recommended spinners for all 5 mm RT probes (P = POM ; K = Kel-F; C = Ceramics):

Sample T [°C]	-150...-80	-80...-50	-50...0	0...80	80...120	120...150
Spinner	C	C	K & C	P, K & C	K & C	C
Recom. VT gas [L/hour]	1200	1000	750	400 (P) 600 (K & C)	450	350
Shim gas [l/min]	20	20	20	0	0...20	20...60
Flush gas [l/min]	5...10	5	5	0	5	5
Chiller	LN2	LN2	BCU II (-50)			

IMPORTANT NOTES (INSTUMENT DAMAGE WILL OCCUR IF NOT FOLLOWED):

1. The temperature of the shim system should always be $-80\text{ °C} < T < 80\text{ °C}$.
2. The temperature of the magnet flange (O-Ring!) has to be $3\text{ °C} < T < 80\text{ °C}$, especially for long term experiments.
3. The shim system and the probe has to be flushed during low temperature experiments (sample temperature $T < 0\text{ °C}$) to prevent icing and condensed water.
4. The new BBFO smart probes on Au400 and Hg400 automatically adjust the Shim/Flush gas based on temperature setting. The A600 needs manually adjustments of Shim/Flush gas to 20/10 from 10/6 LPM (ask NMR staff if you do not know how).
5. You should raise/lower the temperate incrementally by 10 °C for about every 5 minutes
6. Never turn VT gas off.

Table 2: BCU II cooling mode (off, low, medium, and strong)

25 °C to 150 °C	off
25°C down to 0°C	low
0°C to -20°C	medium
-20°C to -50°C	strong

IMSERC User Manual for NMR (V2.0)

Fig 2: Au400/Hg400 VT interface with BCU II (by type command "edte")





Fig 3: A600 VT interface with FTS (by type command "edte")



IMSERC User Manual for NMR (V2.0)

APPENDIX C: 1D SELECTIVE NOESY WITH TOPSPIN

REFERENCE SPECTRUM

Run a 1D Proton spectrum, following the instructions in the short Bruker manual.

SET UP THE 1D SELECTIVE NOESY

The selective pulse regions are set up using the integration tools. Power and duration of the shape pulses are calculated using the hard 90° pulse in the prosol table.

1. Stay in the reference spectrum, on the menu bar, click **Acquire**.

On the **More** button, click the **drop-down** arrow to see more options.




2. In the list, select **Setup Selective 1D Expts.**

The Workflow button bar changes for setting up the 1D selective experiment.

3. On the Workflow button bar, click **1D Selective Experiment Setup**.

4. In the message window, click **Close**.

IMSERC User Manual for NMR (V2.0)



There is no other function to this button then the instruction displayed above.

Expand the spectrum region having peaks you are interested in doing NOE.

5. On the Workflow button bar, click **Define Regions**.  Define Regions

The Define Regions toolbar is displayed:



6. Integrate the peak (multiplet) that you will irradiate to observe the NOE.

If desired, other peaks can be integrated and a separate dataset will be created for each integral saved in the region file.

7. On the toolbar, click **Save/export regions** .

8. In the list, select **Save Regions to 'reg'**.

9. On the toolbar, click **Return do NOT save regions!**. 


10. In the message window, click **No**.



IMSERC User Manual for NMR (V2.0)

11. On the **Create Dataset** button, click the **drop-down** arrow to see more options.

12. In the list, select **Selective gradient NOESY**.




The default parameters are taken from the standard parameter set SELNOGP. The mixing time D8 is dependent on the size of the Molecule and the magnetic strength. It can vary from a large Molecule to a small one from 100 ms to 800 ms. If desired, the Gaus1_180r.1000 pulse can be changed by clicking on the Shape button in the above window. Number of scans (NS) should be determined based on your sample concentration.

13. Enter:

D8 = **0.450**

NS = **32**


14. In the SELNOGP window, click **Accept**.



The new dataset is created and all parameters are automatically set.

IMSERC User Manual for NMR (V2.0)

15. In the sel1d window, click **OK** to start the acquisition.



If you click “cancel”, dataset will be created but not run. You can make further changes and then start acquisition.

DISPLAY 1D NOESY TOGETHER WITH 1D PROTON SPECTRUM

1. Display the selective NOESY spectrum.

2. On the toolbar, click **Multiple display**.  The Multiple display toolbar is displayed:



3. Drag the Reference spectrum (1D proton) into the spectral window.

IMSERC User Manual for NMR (V2.0)

APPENDIX D: PROTON DOSY EXPERIMENT

INTRODUCTION

The **DOSY (Diffusion-Ordered Spectroscopy) experiment** provides accurate, noninvasive, molecular diffusion measurements on biofluids, complex chemical mixtures and multi component solutions. In DOSY spectra, chemical shift is along the detected F2 axis and diffusion coefficient is along the other F1 axis.

Molecules in the solution state move. This translational motion is known as Brownian molecular motion and is often simply called diffusion or self-diffusion. Molecular diffusion depends on a lot of physical parameters like size and shape of the molecule, temperature and viscosity.

Pulsed field gradient NMR spectroscopy can be used to measure translational diffusion. By use of a gradient pulse, molecules can be spatially labeled. After this encoding gradient pulse (δ), molecules move during the diffusion time (Δ). Their new position can be decoded by a second gradient pulse. This encoding/decoding procedure results

$$I(g) = I(0) \exp \left[-(\gamma g \delta)^2 D \left(\Delta - \frac{\delta}{3} \right) \right]$$

in an attenuation of the NMR signal which can be described by the following equation:

Where **I** is the observed intensity, **D** is the diffusion coefficient, γ is the gyro magnetic ratio of the encoded nucleus, **g** is the gradient strength, δ is the length of the gradient pulse, and Δ as mentioned previously is the diffusion time.

The diffusion experiment records a series of 1D ^1H spectra at increasing gradient strengths (**g**) and then fits the signal intensity decay to the above equation to obtain **D**.

Convection within the sample tube, such as, moving liquid columns along the sample axis (primarily due to temperature gradients), can seriously affect diffusion experiments, in particular, at elevated temperatures. Convection currents are caused by small temperature gradients in the sample and result in additional signal decay that can be mistaken for faster diffusion.

DOSY uses three parameters to define the duration of the diffusion: gradient length δ (**P30** in topspin, the total gradient defocusing time), the diffusion gradient level **g** (**GPZ6** in topspin, maximum 95%), and the diffusion delay

IMSERC User Manual for NMR (V2.0)

Δ (**D20** in topspin, 60 ms as default, max determined by the shortest T1 relaxation). In most case, **GPZ6** is the variable parameter to be arrayed for DOSY purpose. Depending on sample, you might need increase **D20** and/or **P30** (max **2ms!**) in order to obtain enough signal attenuation. The purpose of doing this is to get a diffusion decay curve like in the **figure C** below that will give you the best DOSY fitting.

There are two sets of parameter files under user directory:

ledbpgp2s ("longitudinal eddy current delay" LED-bipolar gradients pulse sequence)

dstebpgp3s (double stimulated echo for *convection compensation* and LED using bipolar gradient pulses for diffusion using 3 spoil gradients).

With LED, magnetization is stored along the z-axis during most of the pulse sequence, so T1 relaxation is predominant. Since in macromolecules the T1 relaxation is slower than the T2 relaxation, the LED experiment is better suited to the measurement of D_s of slower diffusing molecules where longer "diffusion delay" is required to detect attenuation of the signal.




Figure: Simulated diffusion decay curves by varying the gradient strength g from 2 to 95% in 16 steps for the same diffusion constant, but with different selection for Δ and δ . They are chosen too small (A), too big (B), and properly (C) to sample data points along the whole decay curve.

EXPERIMENT

The DOSY pulse program used in the following procedure is the Stimulated spin-echo experiment using bipolar gradients and an additional delay just prior to detection for the ring-down of any possible eddy currents (**ledbpgp2s**). The same procedure works for **dstebpgp3s** if you need convection compensation.

IMSERC User Manual for NMR (V2.0)

1. To set up a DOSY experiment, start with recording a normal proton spectrum, followed by optimizing P1, SWH, and O1, if necessary.
2. Type “rpar ledbpgp2s1d_nu all” to retrieve 1D dosy parameters (or “rpar” to select “ledbpgp2s1d_nu”). Update solvent with yours (default is CdCl₃)
3. Check to make sure the P1, SWH, and O1 are same as your proton experiment. The recycle delay D1 should be 1-2 T1 and dummy scan DS should be at least 8. Adjust NS accordingly to give sufficient S/N.
4. Change GPZ6 to 2% and type “zg” to collect data.
5. Use “edc” to create another 1D experiment and change GPZ6 to 75% and type “zg” to collect data
6. Click  (dual display) to compare the 1D data with GPZ6 at 75% to the previous 1D of 2% to check if the nmr signals of interest are attenuated to less than 5-10% of the intensities obtained with GPZ6 at 2%. If you don't get there or already past it, adjust GPZ6 (to 95% or 50%) accordingly to make sure you get there. Write down the GPZ6 value for 2D DOSY setup.
7. If changing GPZ6 alone is not enough to attenuate the signal enough, increase the D20 and/or P30 to achieve the goal.
8. Type “rpar ledbpgp2s_nu all” to retrieve 2D dosy parameters (or “rpar” to select “ledbpgp2s_nu”). Update solvent, P1, SWH, and O1 with the values from your proton experiment
9. Type “dosy” to create the gradient ramp function:
10. Enter first gradient amplitude: 2
Enter final gradient amplitude: 95 (or the value obtained from 1D DOSY)
Enter number of points: 16 (or the number you think appropriate for your sample)
ram type (l/q): l
and finally, Do you want to start acquisition? Select OK to collect 2D DOSY data.

PROCESSING

1. Set the proper window function.
2. Type “**eddosy**”
3. Type “**setdiffparm**” (or click )
4. Type “**xf2**” (or click )
5. If you need phase the spectrum, type “**rser 1**” to read the 1st fid to a new prono and type “**efp**” and “**apk**” to get correct **PHC0** and **PHC1** numbers. Then go back to 2D DOSY dataset and correct the phase values. Remember the phase mode is “**pk**” for direct dimension (F2).

IMSERC User Manual for NMR (V2.0)

6. Type “**dosy2d setup**” (or click )
7. Type “**dosy2d**” (or click ) , you should see the 2D DOSY spectrum with chemical shift along the detected F2 axis and diffusion coefficient along F1 axis.

ADDITIONAL NOTES

Sample preparation: make sure the sample volume is not more than 550 ul.


APPENDIX E: NMR T₁ RELAXATION EXPERIMENT

INTRODUCTION

When an NMR sample sits in the magnet, the applied static magnetic field B_0 will generate the equilibrium magnetization M_0 along +z axis. When a RF pulse is applied to the sample, the net magnetization will be rotated away from +z axis. T₁ relaxation (longitudinal or spin-lattice) is the process by which the net magnetization goes back to its initial maximum value ($M_{z,eq}$) parallel to B_0 .

The inversion-recovery experiment measures T₁ relaxation times of any nucleus. If the net magnetization is placed along the -z axis, it will gradually return to its equilibrium position along the +z axis at a rate governed by T₁. The equation governing this behavior as a function of the time t after its displacement is:

$$M_z(t) = M_{z,eq} \left(1 - 2e^{-t/T_1} \right)$$



The basic pulse sequence consists of an 180° pulse that inverts the magnetization to the -z axis. During the following delay, relaxation along the longitudinal plane takes place. Magnetization comes back to the original equilibrium z-magnetization. A 90° pulse creates transverse magnetization. The experiment is repeated for a series of delay values taken from a variable delay list. A 1D spectrum is obtained for each value of vd and stored in a pseudo 2D

dataset. The longer the recycle delay ($d1$) is, the more precise the T₁ measurement is. Ideally $d1$ should be set to $5 \cdot T_1$. A rough estimation of the T₁ value can be calculated from the null-point value by using $T_1 = t_{null} / \ln 2$.

SETTING UP PROTON T₁ EXPERIMENT

1. To set up a T₁ experiment, start with recording a normal proton spectrum to adjust the spectral sweep width **SWH**, acquisition time **aq** and other parameter if necessary.
2. Create new dataset and load “**Proton_T1**” parameter set. Update the parameters with the ones you obtained from last step. The recycle delay **D1** should be $\sim 2-5 \cdot T_1$. Adjust **NS** accordingly to give sufficient S/N (**fig 1**).
3. Edit the “t1delay” by clicking on  at VDLIST line in **fig 1**. **Fig 2** is a good starting list.

IMSERC User Manual for NMR (V2.0)

4. Change the “TD” value for F1 dimension to the number in your VDLIST (**fig 3**)
5. Collect the pseudo 2D T1 dataset

Fig 1. ACQUPARS display in “pulse program parameters” view

Parameter	Value	Description
PULPROG	t1ir	Pulse program for acquisition
TD	32786	Time domain size
SWH [Hz, ppm]	8196.72 / 13.6657	Sweep width
AQ [sec]	1.9999460	Acquisition time
RG	64	Receiver gain
DW [µsec]	61.000	Dwell time
DE [µsec]	20.00	Pre-scan-delay
D1 [sec]	5.000000000	Relaxation delay: 1-5 * T1
d11 [sec]	0.0299999993	Delay for disk I/O [30 msec]
DS	0	Number of dummy scans
NS	2	Scans to execute
VDLIST	t1delay	Variable delay list
vd [sec]	5.00000000	vd[10]= { 5.000000 sec 0.001000 sec... }

Fig 2. An example of t1delay list with 8 delays

Line	Delay Value
1	0.050
2	0.100
3	0.250
4	0.500
5	0.800
6	1.5
7	3
8	5

IMSERC User Manual for NMR (V2.0)

Fig 3. ACQPARS display in “all acquisition parameters” view

The screenshot displays the ACQPARS window in the "all acquisition parameters" view. The window title is "Probe: CP QCI 600S3 H/F-C/N-D-05 Z". The interface is organized into several sections:

- Navigation Bar:** Includes tabs for SPECTRUM, PROCPARS, ACQPARS (selected), TITLE, PULSEPROG, PEAKS, INTEGRALS, SAMPLE, STRUCTURE, PLOT, FID, and ACQU.
- Left Panel:** A vertical list of parameter categories: Experiment, Width, Receiver, Nucleus, Durations, Power, Program, Probe, Lists, NUS, Wobble, Lock, Automation, Miscellaneous, User, and Routing.
- Main Panel:** A table of parameters with columns for F2, F1, and Frequency axis. The parameters are:
 - Experiment: t1ir (Current pulse program)
 - AQ_mod: DQD (Acquisition mode)
 - FnTYPE: traditional(planes) (nD acquisition mode for 3D etc.)
 - FnMODE: QF (Acquisition mode for 2D, 3D etc.)
 - TD: 32786 (Size of fid)
 - DS: 0 (Number of dummy scans)
 - NS: 2 (Number of scans)
 - TD0: 1 (Loop count for 'td0')
 - TDav: 0 (Average loop counter for nD experiments)

IMSERC User Manual for NMR (V2.0)

PROCESSING

1. Process and adjust phase for the dataset. Use **rser n** (n is the number of total delays) to read out the last fid. Process and phase correct it. On the Adjust Phase toolbar, click **Save for spectrum**.




2. Go back to pseudo 2D T1 dataset by closing the 1D window
3. At the command prompt, type **xf2** to process only the F2 axis. Type **abs2** to baseline correct the rows.
4. On the menu bar, click **Applications**.
5. On the **Dynamics** button, click the drop-down arrow to see more options and in the list, select **T1/T2 Module**.




6. The flow buttons change to determine the T1 / T2 relaxation times. While executing the steps below, message windows will be displayed. Please read each message thoroughly and follow the instructions. On the Workflow button bar, click **Fid**


7. In the Extract a row from 2d data window, click **Spectrum**



8. Enter Slice Number = **n** (the last one).




9. On the Workflow button bar, click **Peaks/Ranges**.



10. In the Define Peaks and/or Integrals window, click **Manual Integration**.

IMSERC User Manual for NMR (V2.0)




11. Define the regions by drawing an integral over the peaks of interest, On the Integration toolbar, click **Save/export integration regions**



12. In the list, select **Export Region To Relaxation Module**.



13. In the Prepare relaxation data window, click **OK**



14. On the Workflow button bar, select **Relaxation**.




15. By default, the selected areas are peak-picked, and the first peak is displayed in the Relaxation window.

16. On the Workflow button bar, select **Fitting**



17. In the message window, click **Close**.



IMSERC User Manual for NMR (V2.0)

APPENDIX F: KINETICS / REACTION MONITORING

INTRODUCTION

As an intrinsically quantitative analytical technique, NMR spectroscopy can be used to measure concentrations of different components during chemical reactions for period of a few minutes to multiple days. One can either take a series of 1D spectra or acquire the data in pseudo 2D mode over the period of reaction. The best way to analyze the resulting data is to use MNOVA.

Before starting, you should have good estimates of the timeframe for your reaction and T1 relaxation time for your sample. Setup the **d1** and **aq** accordingly to make sure the relaxation time (**d1+aq**) is at least $5 \cdot T_1$. In general, the aliphatic protons on small molecules have T1 around 2 s and aromatic proton at about 4 s. You should run T1 experiment to get the numbers since it is very sample/solvent dependent.

1D method is a very versatile. Setup the 1D you will repeat and run **multi_zgvd**. It can be either a single pulse or single pulse with decoupling, i.e. F19 or F19CPD. The drawback of this method is that the timing between spectra can be off by a few seconds especially when your kinetics are fast and take less than an hour. The pseudo 2D, on the other hand, gives perfect timing for each fid, but it can be used only for single pulse experiment.

Once you decide how often you take a spectrum, you need strike a balance between the time resolution of the kinetic measurement and the amount of time needed to obtain sufficiently good signal-to-noise for each experiment. Limit the number of scans (ns) to be as small as necessary for adequate signal-to-noise to improve time resolution.

Before starting your reaction, please setup the experiment you want to repeat with a test sample with conditions similar to your real one. Do the locking, tuning, and shimming. Find the appropriate number of scans (ns) for adequate signal-to-noise. If resolving peaks is not a concern, you do not have to do lock/tune/shimming after putting the real sample in. Simply start acquiring the data, especially for F19.

SEPARATE 1D SPECTRA W/ MULTI_ZGVD

This works for any nucleus, including proton.

1. Assume you already determined how often you want to run your recurring 1D experiment, let's call it **D20**, delay between start of different 1D spectra

IMSERC User Manual for NMR (V2.0)

2. Start with a normal 1D spectrum to adjust the spectral sweep width **SWH**, acquisition time **aq**, offset **O1p**, number of scans **NS**, and other parameter obtain sufficiently good signal-to-noise if necessary. Type **expt** to calculate how much time it takes. Let's call it "**T_{expt}**". The delay between the end of one fid and start of next one equals **D20-T_{expt}**. Let's call it **D_{fix}**.
3. Create a new dataset with exactly same parameters from step 1. Start your reaction and load your sample to NMR instrument as fast as you can. Since you have already done locking/tuning/shimming on a test sample with similar conditions, you have following options:
 - a) Do a topshim session first if your kinetics takes hours to finish.
 - b) Skip the topshim if your reaction is really fast
4. Run **multi_zgvd**, when asked for a fixed or variable delay, answer with the default (fixed delay), then give the **D_{fix}** as the input for next question. For the question of "Enter number of experiments", give the numbers of experiment you want to run.
5. During the run, you can use **multiple display** to check peak intensity changes to evaluate if you reaction finishes or not.

PSEUDO 2D MODE PROCEDURE

The following procedure can be used for any nucleus.

1. Following step 1 and 2 of previous section to optimize the 1D experiment you want to repeat.
2. Create a new dataset and load the parameter set "**kx_zg2d_nu**".
3. Input the **D20** (delay between start of different 1D spectra) as shown in **Fig 1** and **TD** on F1 dimension (how many 1D spectra you want to acquire) as shown in **Fig 2**.
4. Start your reaction and load your sample to NMR instrument as fast as you can. Since you have already done locking/tuning/shimming on a test sample with similar conditions, you have following options:
 - a) Do a topshim session first if your kinetics takes hours to finish.
 - b) Skip the topshim if your reaction is really fast
5. Start your experiment by typing **zg** or click on "**run**".
6. During the run, you can use **rser** to check each individual fid as long as it is finished. For example, "**rser 1 10**" will write the 1st fid to experiment number 10; "**rser 20 11**" will write the 20th fid to experiment number 11. Then you can use **multiple display** to stack or superimpose them.

IMSERC User Manual for NMR (V2.0)

Fig 1. ACQUPARS display in “pulse program parameters” view

7 test 51 1 /home/walkon/data/Zhang/yzh933

Spectrum ProcPars **AcquPars** Title PulseProg Peaks Integrals Sample Structure Plot Fid Acqu

Probe: PA BBO 400S1 BBF-H-D-05 Z SP N

General Channel f1

General

PULPROG lx_zg2d_nu Pulse program for acquisition

TD 25606 Time domain size

SWH [Hz, ppm] 6393.86 15.9958 Sweep width

AQ [sec] 2.0023892 Acquisition time

RG 18 Receiver gain

DW [µsec] 78.200 Dwell time

DE [µsec] 6.50 Pre-scan-delay

D1 [sec] 10.000000000 Relaxation delay; 1-5 * T1

D20 [sec] 276.000000000 Delay between start of different 1D spectra

D21 [sec] 0 Shift delay for the first increment

DELTA [sec] 179.95080566 DELTA=d20-((d1+p0+de+aq)*(ns+ds))-30m

DS 0 Number of dummy scans

NS 8 1 * n, total number of scans: NS * TD0

ZGOPTNS Options for zg

Channel f1

SFO1 [MHz] 399.7218787 Frequency of ch. 1

O1 [Hz, ppm] 1878.68 4.700 Frequency of ch. 1

NUC1 1H Edit... Nucleus for channel 1

CNST18 30.0000000 Flip angle in degree

p0 [µsec] 3.33 For any flip angle

P1 [µsec] 10.000 F1 channel - 90 degree high power pulse

PLW1 [W, dB] 15.162 -11.81 F1 channel - power level for pulse (default)

Fig 2. ACQUPARS display in “all acquisition parameters” view

7 test 51 1 /home/walkon/data/Zhang/yzh933

Spectrum ProcPars **AcquPars** Title PulseProg Peaks Integrals Sample Structure Plot Fid Acqu

Probe: PA BBO 400S1 BBF-H-D-05 Z SP N

Experiment Width Receiver Nucleus Durations Power Program Probe Lists NUS Wobble Lock Automation Miscellaneous User Routing

Experiment

PULPROG lx_zg2d_nu Current pulse program

AQ_mod DQD Acquisition mode

FnTYPE traditional(planes) nD acquisition mode for 3D etc.

FnMODE QF Acquisition mode for 2D, 3D etc.

TD 25606 16 Size of fid

DS 0 Number of dummy scans

NS 8 Number of scans

TD0 1 Loop count for 'td0'

TDav 0 Average loop counter for nD experiments

Width

SW [ppm] 15.9958 10.0000 Spectral width

SWH [Hz] 6393.862 3997.219 Spectral width

IN_F [µsec] 250.17 Increment for delay

AQ [sec] 2.0023892 0.0020014 Acquisition time

FIDRES [Hz] 0.499403 499.652344 Fid resolution

FW [Hz] 4032000.000 Filter width

IMSERC User Manual for NMR (V2.0)

APPENDIX G: EVANS METHOD

INTRODUCTION

Evans Method was established in 1959 for magnetic susceptibility measurement [Ref 1]. It was developed by many other researchers in the past half century.

This Evans Method manual will help NMR users understand paramagnetic susceptibility measurement by using NMR. Users can simply follow the procedures to conduct their experiments with little or even no NMR staff assistance

Now, the most common equation is (1) below

$$\chi_{mass} = 3\Delta f/4\pi f m + \chi_0 + \chi_0(d_0 - d_s)/m \quad (1) \text{ [Ref 2]}$$

- χ_{mass} is mass Susceptibility in cm^3g^{-1}
- Δf is obs freq diff in Hz
- f is spectrometer freq, e.g. 399.732×10^6 Hz on Hg400
- m is mass of paramagnetic substance in $\text{g}\cdot\text{cm}^{-3}$ i.e. concentration in g/mL
- χ_0 is mass susceptibility of solvent in cm^3g^{-1}
- d_0 is density of solvent in $\text{g}\cdot\text{cm}^{-3}$
- d_s is density of solution in $\text{g}\cdot\text{cm}^{-3}$

There are some other related concepts: [Ref 3]

volume Susceptibility (χ_v) $\chi_v = \rho\chi_{mass}$ here ρ is the density in $\text{kg}\cdot\text{cm}^{-3}$ or $\text{g}\cdot\text{cm}^{-3}$ and

molar Susceptibility (χ_{mol}) $\chi_{mol} = M\chi_{mass}$ here M is molecular weight

Also, understand Conversion of SI and cgs units: $\chi^{SI} = 4\pi\chi^{cgs}$

In this manual, we use cgs unit through-out. For instance, χ_{mol} in $\text{cm}^3\text{mol}^{-1}$


PREPARATION

Prerequisite: users have done the basic NMR training

Spectrometer: NMR-Au400 recommended

NMR tube: coaxial inserts (see the picture below) or capillary inserts [Ref 4]

IMSERC User Manual for NMR (V2.0)



NMR Sample Prep: (using FeSO₄, t-butanol and D₂O as example [Ref 5])

- Solution A: dissolve 7.0 mg FeSO₄·7H₂O in 0.5 mL D₂O, add 30 μ L t-butanol, and adjust total volume to 1 mL with D₂O
- Solution B: 30 μ L t-butanol + 970 μ L D₂O
- Transfer A into the coaxial insert and B into a 5mm tube

EXPERIMENT SETUP (QUICK PROCEDURES)

1. Turn on spectrometer computer monitor by login NUcore with netid/password
2. Create a new data set and Set up a regular 1D 1H experiment
3. Insert NMR sample into magnet with SX
4. Do locking, tuning and shimming as usual
5. Run a regular 1D H1 (will see two methyl peaks caused by paramagnetic reagent)

DATA ANALYSIS

1. measure distance between two methyl peaks (Δf 510 Hz)
2. use a simplified equation $\chi_{mass} = 3\Delta f/4\pi f m + \chi_0$ (2)

Here, $f = 399.732$ on Hg400

$m = 0.0072$ g (actual weight) in 1.00 mL (solution A)

χ_0 is approximately equal to water mass susceptibility, -0.72×10^{-6} cm³/g

3. molar susceptibility $\chi_{mol} = M\chi_{mass}$
Here, molar mass M is 278 for FeSO₄, and $\chi_{mol} \sim 11200 \times 10^{-6}$ [Ref. 6]
4. Actual experimental result is $\chi_{mass} = 41.58 \times 10^{-6}$ cm³/g and $\chi_{mol} = 11600 \times 10^{-6}$ cm³/mol
5. If the χ_{mol} value is small, you may consider the diamagnetic susceptibility correction term

$$\chi_{mol\text{-para}} = \chi_{mol} + \chi_{mol\text{-dia}} \text{ [ref 7]}$$

IMSERC User Manual for NMR (V2.0)

REFERENCES

1. D.F. Evans, J. Chem. Soc. 2003 (1959)
2. <http://pubs.acs.org/doi/pdfplus/10.1021/ed069p62.1>
3. http://en.wikipedia.org/wiki/Magnetic_susceptibility
4. <http://www.wilmad-labglass.com/Support/NMR-and-EPR-Technical-Reports/NMR-007--Coaxial-Inserts-in-NMR-Studies/>
5. 200 and More NMR Experiments. Stefan Berger and Siegmund Braun (2004)
6. http://www-d0.fnal.gov/hardware/cal/lvps_info/engineering/elementmagn.pdf
7. <http://pubs.acs.org/doi/pdf/10.1021/ed085p532>

IMSERC User Manual for NMR (V2.0)

REVISIONS

V1.0 2020/2/2	<ul style="list-style-type: none">• Initial release.
V1.1 2020/2/10	<ul style="list-style-type: none">• Added 'Troubleshooting' section
V2.0 2020/4/28	<ul style="list-style-type: none">• Added additional specialty NMR experiments
