

IMSERC PULSE Q-band Manual

CONTENTS

Introduction.....	2
Safety.....	2
Data Management.....	2
Software	3
Materials & methods.....	3
Experiment Setup	3
Step 1. Sample preparation.....	3
Step 2. Turn on the spectrometer and open the software	4
Step 3. Insert sample.....	5
Step 4. Setting temperature.....	7
Step 5. Connect the spectrometer and create new experiment.....	7
Step 6. Tuning the spectrometer.....	8
Step 7. Checking Defense pulse.....	10
Step 8. Turn on TWT.....	14
Step 9. Optimize signal	15
Step 10. Optimize signal with a nutation experiment.....	17
Step 11. Optimize signal with the bandwidth experiment.....	18
Step 12. Collecting 1D EPR spectra.....	19
Step 6. Saving the data	21
Ending Work	21
Troubleshooting	22
Revisions.....	22
Appendix A: How to run with pulse sequence EPR	23
Appendix B: T_1 and T_2 measurement.....	25
Appendix C: ELDOR(DEER) experiment	27

IMSERC Pulse Q-band Manual (V1.0)

INTRODUCTION

The Bruker pulse Q-band EPR spectrometer provides high-field pulsed EPR capability for advanced characterization of paramagnetic systems. Operating in the Q-band frequency range, the instrument offers improved spectral resolution compared with X-band EPR, particularly for resolving g-anisotropy and distinguishing overlapping paramagnetic species. The system is well suited for detailed studies of transition-metal complexes, radicals, and spin-labeled biomolecules. With cryogenic temperature capability, the Bruker pulse Q-band EPR enables low-temperature pulsed EPR experiments, including echo-detected field-swept EPR, relaxation measurements, and DEER for determining nanoscale distance distributions.

SAFETY

All users of IMSERC must review the general safety policies at <http://imserc.northwestern.edu/about-policies.html>. Plus, users should know the points below:

- Users with pacemaker are not allowed
- Emergency Exits
- Location of shower and eyewash
- No food and drink
- Dressing code

DATA MANAGEMENT

IMSERC Pulse Q-band Manual (V1.0)

The personal data folder is created during training. The data were under the user's personal folder, which must be located under the supervisor's group folder.

EPR data on this instrument are copied on 'imsercdata.northwestern.edu' under 'data2/pcm/Qband' within a few seconds. Please follow instructions at <http://imserc.northwestern.edu/about-general-faq.html#data> for details about data access.

SOFTWARE

The EPR experiment on this instrument is performed with Bruker Pulse Q band.

MATERIALS & METHODS

The default measurement mode of pulse Q band with the 1.6 mm Bridge 12 resonator. If user want to run an experiment in a different resonator, they will contact the staff. The staff will also need to tune back to the default resonator after the experiment. (these will be charged via staff time)

The default working condition of INSTRUMENT-NAME is as follows:

1. Computer screen is by default deactivated. Users must start the reservation through NUCore in order to be able to turn on the computer screen.
2. Leave the acquisition software open when user complete with the measurement.

EXPERIMENT SETUP

STEP 1. SAMPLE PREPARATION

Sample Tubes:



IMSERC Pulse Q-band Manual (V1.0)

a) Standard EPR tube size = 1.6 mm OD, however, sizes do vary (3mm, 5mm) based on different resonators. 1.2 mm OD capillary tube can be used a second tube to avoid sample explosion

b) EPR tubes are made out of quartz whereas NMR tubes are made out of Pyrex, which contains some iron that may or may not, depending on the sample, contribute to the EPR signal.

Tubes can be purchased by various vendors. A couple of options are,

Wilmad Glass: <https://sp-wilmadlabglass.com/product/q-band-cfq-epr-tube-100mm-l-one-end-sealed-10-pk/>

Aqueous Samples (also true for highly polar and conducting samples): The water in aqueous samples absorbs the electric part of the microwaves and begins to heat (like your microwave at home), and therefore, the sample size needs to be limited, specifically the surface area of the sample size. Current, IMSERC use capillary tubes to adsorb the liquid sample and wax to seal each open side. Then input the capillary tube with sample into the 1.6 mm quartz tube. The concentration of liquid sample is about 100 – 200 μ M.

The liquid sample should be also pre-freeze in capillary tube to avoid the explode.

Powder samples: the capillary tube was sealed firstly, and power was packed into the tube.

STEP 2. TURN ON THE SPECTROMETER AND OPEN THE SOFTWARE

1) Turn on the following (if not already on):

Console (if it is on, then skip)

Power Supply of the magnet (AC, reset, then DC)

IMSERC Pulse Q-band Manual (V1.0)



2) turn on the reservation on Nucore. Open the Acquisition software “XEPR” located on the desktop.

STEP 3. INSERT SAMPLE


1) Added around 1 cm sample to the 1.6 mm OD tube. Then assemble the 1.6 mm OD tube with the sample holder. We only need maximum 1 cm sample. The distance from sample holder to the center of sample is about 1 inch. Connect the sample holder to the sample rod.

IMSERC Pulse Q-band Manual (V1.0)



2) Wipe the outside of the tube with acetone or isopropanol.

3) Loosen the plug cap in the top of resonator and remove the whole plug. Quickly insert the sample rod with sample into the resonator. The whole stuff will go straight down to the end and then tighten the metal screw nut to seal the rod to the resonator.



IMSERC Pulse Q-band Manual (V1.0)

STEP 4. SETTING TEMPERATURE

The temperature should be optimized for each sample type, but 40K should be a good starting approximation. – control, set temperature, auto PID

STEP 5. CONNECT THE SPECTROMETER AND CREATE NEW EXPERIMENT

1) Open XEPR.


2) Go to connect to the spectrometer. The led light on the bridge should be on.

Where: “Acquisition” → “Connect to spectrometer” → “local host” → “ok”


3) Create a new experiment and set parameters to hardware.

Where: “Create Experiment” → “pulsed” → “Advanced” “create”

Click “send parameters to hardware”



IMSERC Pulse Q-band Manual (V1.0)



STEP 6. TUNING THE SPECTROMETER

Go to the “Bridge Tuning”, the attenuation should be at 60 dB.

In the tune mode adjust the attenuation to 15dB and find the dip. The true dip should respond to the sample movement and the iris screw rotation (changing the resonator coupling).

Then change from “tune” to “operate” mode. At this moment, the frequency will be set to the right value.


Make sure the resonator is over coupled, this will give us the better performance due to the large bandwidth.

When the resonator is critically coupled, the bandwidth is about 80 mHz



When the resonator is over-coupled, the bandwidth is about 400 mHz.

Overcoupling the resonator also prevent receiver damage due to intense reflections.

IMSERC Pulse Q-band Manual (V1.0)




Critical coupled condition:



The iris should be turn up until the dip almostly disappears, then you achieve the over-coupled condition.

If the resonator is critical-coupled and you will find a significant ring-down effect when you check the echo signal later. Such as this. If you see this, please change the iris to over-coupled.

IMSERC Pulse Q-band Manual (V1.0)




STEP 7. CHECKING DEFENSE PULSE

- 1) Switch to Pulse mode.

IMSERC Pulse Q-band Manual (V1.0)


Where: "FT Bridge" → "Bridge Configuration" → from "CW mode" to "pulse"




2) Check for the defense pulses.

Open "Parameters button" → "Patterns"

Under channel "Acquisition Trigger", program a 20 ns (length, first column) acquisition trigger, position as "0". Change integration time to 1.




IMSERC Pulse Q-band Manual (V1.0)




Under channel "AWG1", program echo detection pulse.

	Pulse 1	Pulse 2
Position	0 ns	200 ns
Length	6 ns	12 ns

IMSERC Pulse Q-band Manual (V1.0)



Open “FT parameters” → “Microwave”, Adjust video bandwidth to 20 MHz by arrow button and video gain to ~ 30 db. Keep the attenuation at 60 dB.




Then clicked the “send parameters to hardware” button in the left bottom corner of XEPR.

Click “Start” under “Patterns”. Open “SpecJet” and click “Run”. Set Number of averages in SpecJet to 100. Now we should see the defense pulses, we may need to adjust the number of points, such

IMSERC Pulse Q-band Manual (V1.0)

as 2048 for that. Make sure you clicked the “send parameters to hardware” button to enable the pulse operation.



STEP 8. TURN ON TWT

If you see the defense pulses then you can turn on the TWT. To turn on the TWT, please press the on/off, wait until the TWT shows “standby”. When TWT is “standby”, press “operate” to turn on the HV (-17.9 V). But if the defense pulse is not visible then do not turn on TWT.

IMSERC Pulse Q-band Manual (V1.0)




STEP 9. OPTIMIZE SIGNAL

- 1) Decrease Attenuation "FT EPR Parameters" → "Microwave" → "attenuation" OR the "Receiver unit" in "FT bridge"

IMSERC Pulse Q-band Manual (V1.0)


- a. You should see an echo start to form. You may have to adjust the video gain to decrease the signal as you go. The signal will initially become stronger when you lower the attenuator then usually goes through a maximum and starts to decrease., Find the attenuation value for the maximum signal.
- b. If you do not find an echo, you possible miss the correct EPR signal frequency and the magnetic field. So, you need to adjust the frequency and the corresponding magnetic field to find the echo signal. If you do not see an echo, you possible miss the correct value of the magnetic field for your ESR signal. In this case, you need to adjust the magnetic field value to find the echo signal.

For a usual nitroxide sample with a g factor of ~ 2.006 at the frequency of $\nu=34$ GHz, the magnetic field value is ~ 12110 Gauss. The magnetic field position is related to the ESR frequency as $B=h\nu/g\beta$, where $h=6.626 \times 10^{-27}$ erg/s is the Planck and $\beta=9.274 \times 10^{-21}$ erg/G is the Bohr magneton.



- 2) When the echo signal is found, change the position of the acquisition trigger to the beginning of the echo, ~ 1140 , to skip the defence pulses. Reduce the no of points to 256 and you will observe the signal only.

IMSERC Pulse Q-band Manual (V1.0)



- 3) Phase the echo using “FT Bridge” → “Receiver Unit” → “Signal Phase”, make all the signal go to the “green curve”, and make the “yellow curve” small, almost flat
Adjust the field position to the maximum of the signal.

Adjust the video gain so that the echo does not exceed the specjet windows.

Adjust the attenuation to the maximal signal, usually between 1 and 7 dB depending on the sample size, sample properties, resonator coupling and temperature.

Use arrow button with ctrl and change the gain/attenuation by 10 directly.

STEP 10. OPTIMIZE SIGNAL WITH A NUTATION EXPERIMENT

- Under “FT Bridge” → “Acquisition” → “Acquisition Mode”, change to “Read Transient” and then click the “play” button.
- Determine the width of the echo to integrate over and put this number into the “Length” of the acquisition trigger.
- Under “Acquisition Trigger”

IMSERC Pulse Q-band Manual (V1.0)

- i. Add 4000 to "Position".
 - ii. Set "Pos. Disp." To 2 ns.
- d. Under "AWG1"
 - iii. Click "Edit" and "Insert Column"
 - iv. Add 4000 ns to the position of pulse 2 and 3.
 - v. Set pulse 1 "Length Increment" to 2 ns.
 - vi. Set "Position Displacement" of pulse 2 and 3 to 2 ns.

The screenshot shows the 'FT EPR Parameters' window with the 'Patterns' tab selected. Under 'PULSE PATTERNS', 'Channel Selection' is set to 'AWG1'. 'Shot Rep. Time [us]' is 1020.00 and 'Shots Per Point' is 100. An 'Edit' button is visible. Below it is a table with 5 columns and 5 rows of parameters. The 'Pos. Disp. [ns]' row has a value of 2 in the 3rd column, which is highlighted with a box and a small up/down arrow. To the right of the table are 'Start' and 'Stop' buttons.

	1	2	3	4	5
Position [ns]	0	4000	4200	0	0
Length [ns]	0	7	14	0	0
Pos. Disp. [ns]	0	2	2	0	0
Length Inc. [ns]	2	0	0	0	0
Frq. Start [MHz]	0	0	0	0	0

- e. Under "FT Bridge" → "Patterns" set "Shots per point" to 10
 - f. Under "FT Bridge" → "Acquisition" set "X-axis Size" to 256
 - g. Change X-Axis Quantity to "Time"
 - h. Change Acquisition mode to "Run from Tables"
 - i. Click "start"

STEP 11. OPTIMIZE SIGNAL WITH THE BANDWIDTH EXPERIMENT


Run bandwidth experiment.

- a. Set MW frequency to 34 GHz and magnetic field at 12080 Gauss
- b. Load Prodel program
 - vii. "Processing" → "ProDel" → "Load Programs"
 - viii. /home/xuser/xepFiles/ProDel/sharedProDel/Standards

IMSERC Pulse Q-band Manual (V1.0)

- ix. Switch to show "FileNames"
- x. Load "BndWidth_Nutation_Qstandalone.ppp"
- c. Run the program by selecting it under "ProDel".
 - xi. Number of points: 40
 - xii. Step Size: 10
 - xiii. Tolerance: 1

After the bandwidth experiment, you will find the best frequency and magnetic field at maximum amplitude, which will give the highest efficiency of RF pulse (shortest pi pulse). However, make sure check the nutation curve shape at the good amplitude.




STEP 12. COLLECTING 1D EPR SPECTRA

- (1) Set the magnetic field and frequency to the number found in the bandwidth experiment.
- (2) Under "FT parameters" → "Acquisition" set Acquisition mode to "ran from PulseSpel" (load program: IMSERC_Qband.exp, load parameter from IMSERC_Qband.def)

IMSERC Pulse Q-band Manual (V1.0)

- (3) In the “experiment” chose “2P ESE set up”, phase cycle “ 2 step”. Make sure the set p0 and p1, attenuation and gain from previous steps.
- (4) In pulseSPEL variable, type d0 and set it to 0.
- (5) Type d1 and set it to 180 ns.
- (6) Type srt, set it to 1000*1.02 ns, compile and validate the variables and program.
- (7) Click the “run” in the left bottom corner and record the echo signal. Measure the time position of the point at the beginning of the echo signal which is located at approximately $\frac{1}{4}$ of the hight.
- (8) After that, change the d0 to adjust the beginning of echo to this point.
- (9) Measure the echo signal width at approximately $\frac{1}{4}$ height.
- (10) Type pg in pulseSPEL variable and set this number for integration of echo signal, for example 74 for the nitroxide echo.
- (11) Switch pulse sequence from “2p ESE set up” to “ 2P ESE field sweep”
- (12) Under “FT parameters” $\rightarrow$ “Acquisition” set x axis quantity “Magnetic field”
- (13) Under “FT parameters” check the filed center and sweep width.
- (14) Set the correct dim number in program.
- (15) Click the “run” in the “left bottom corner”. You should get a 1D ESR absorption spectrum
- (16) Click “save” under “file” to save the spectra data.

IMSERC Pulse Q-band Manual (V1.0)



STEP 6. SAVING THE DATA

As a File

File => Save as... => My Computer => D Disk => User's name folder

As an ASCII file

File => Export... => Destination (notice the file type is ASCII format)

Open this file in Excel and it is delimited with tabs; format as necessary

ENDING WORK

- 1) Change the attenuation to 60 dB

IMSERC Pulse Q-band Manual (V1.0)

- 2) Set the bandwidth back to 200 MHz
- 3) Before you turn off the TWT, we need to run the "TWT maintenance". Change the pulse sequence to "twt maintenance".
- 4) Set n=100
- 5) Click "play". When this happens, the TWT is on the operate mode.
- 6) After 20 mins, stop the "TWT maintenance".
- 7) change the TWT to standby by pressing the standby button, wait until it achieves standby condition
- 8) change the mode from "pulse" to "CW" (FT bridge)
- 9) In the tuning bridge, change to "tune" then change it to "standby"
- 10) Turn off the magnetic field off by pressing DC then AC
- 11) Disconnect the spectrometer (acquisition tab in software)
- 12) Close the XEPR software
- 13) Change the setting of Temperature back to 45 K.
- 14) Move your sample (take out from the resonator). Put the black holder back to the resonator and tighten it.
- 15) Keep desk top clean
- 16) Logout from NUcore

TROUBLESHOOTING

For any problems in pulse Q-band, please contact Dr. Jinlei Cui and Dr. Chris Malliakas

REVISIONS

V 1.0 2025/02/16	Release of original version of IMSERC Pulse Q-band V1.0
---------------------	---

IMSERC Pulse Q-band Manual (V1.0)

APPENDIX A: HOW TO RUN WITH PULSE SEQUENCE EPR

For many experiments, such as T1, T2, ELDOR and ESEEM, it is more convenient to run from the pulse sequence.

For example, you can load the parameters and pulse sequence by the following methods for DEER.

- (1) Click “FT Bridge” → “Acquisition” → “PulseSPEL”
- (2) Load Program
 - xiv. /home/xuser/xexprFiles/PulseSPEL/IMSERC/IMSERC_Qband.exp
- (3) Load Variable Definitions
 - xv. /home/xuser/xexprFiles/PulseSPEL/IMSERC/IMSERC_Qband.def
- (4) Compile variable definitions and compile/validate the program.

For 1D echo field sweep, T1 and T2 relaxation.





- Load Program
 - i. /home/xuser/xexprFiles/PulseSPEL/IMSERC/relaxation.exp
- Load Variable Definitions
 - ii. /home/xuser/xexprFiles/PulseSPEL/IMSERC/relaxation.def

For ESEEM.

- Load Program
 - i. /home/xuser/xexprFiles/PulseSPEL/IMSERC/3p_eseem.exp
- Load Variable Definitions
 - ii. /home/xuser/xexprFiles/PulseSPEL/IMSERC/3p_eseem.def

The difference between relaxation.exp and DEER.exp is the former one uses the spg channel and the latter one uses the awg channel.

IMSERC Pulse Q-band Manual (V1.0)



IMSERC Pulse Q-band Manual (V1.0)

APPENDIX B: T_1 AND T_2 MEASUREMENT

For T_1 and T_2 measurement, please first follow the previous steps to acquire the 1D EPR spectra. Then follow Appendix A to load the pulse sequence and parameters. Then please take the following steps to run the T_1 and T_2 measurement.

T_2 measurement:

- Under “FT parameters” → “Acquisition” set Acquisition mode to “ran from PulseSpel”
- In the “experiment” chose “2P ESE set up”, phase cycle “2 step”. Make sure the set p0 and p1, attenuation and gain from previous steps.
- In pulseSPEL variable, type d0 and set it to 0.
- Type d1 and set it to 180 ns.
- Type srt, set it to 1000×1.02 ns
- Click the “run” in the left bottom corner and record the echo signal. Measure the time position of the point at the beginning of the echo signal which is located at approximately $\frac{1}{4}$ of the height.
- After that, change the d0 to adjust the beginning of echo to this point.
- Measure the echo signal width at approximately $\frac{1}{4}$ height.
- Type pg in pulseSPEL variable and set this number for integration of echo signal, for example 74 for the nitroxide echo.
- Switch pulse sequence from “2p ESE set up” to “2P ESE decay”
- Type d30 and set it to 4 ns.
- Change the dim2 to correct number for the echo decay. It should allow for seeing the whole echo decay to nearly zero.
- Make sure P0 and P1 are correct.
- Compile variable definitions and validate the program.
- Change the phase cycle to “2 steps”
- Click the “run” in the “left bottom corner”. You should get an ESR echo decay curve.

T_1 measurement:

- Load the program IMSERC/relaxation.exp
- Do the same steps from step a to step h in T_2 measurement.
- Switch pulse sequence from “2p ESE set up” to “saturation recovery”
- Type d30 and set it to 2000 ns. d19 to 8000 ns and d2 to 750 ns.
- Change the dim5 to correct number for the saturation. It should allow for seeing the whole saturation recovery.

IMSERC Pulse Q-band Manual (V1.0)

- f. Make sure P0 and P1 are correct.
- g. Set $srt > 2 * d30 * dim5$
- h. Compile variable definitions and validate the program.
- i. Change the phase cycle to “2 steps”
- j. Click the “run” in the “left bottom corner”. You should get an saturation recovery curve.

IMSERC Pulse Q-band Manual (V1.0)

APPENDIX C: ELDOR(DEER) EXPERIMENT

For ELDOR measurement, please first follow steps in main manuscript to run the 1D echo. However, you need to use the awg channel instead of the spg channel.

Then following the below steps to run the ELDOR optimization and ELDOR measurement.

Click “FT Bridge” → “Acquisition” → “PulseSPEL”

Load Program

xvi. /home/xuser/xexprFiles/PulseSPEL/IMSERC/IMSERC_Qband.exp

Load Variable Definitions

xvii. /home/xuser/xexprFiles/PulseSPEL/IMSERC/IMSERC_Qband.def

Compile variable definitions and compile/validate the program.

1. T2 measurement

- a. Under “FT parameters” → “Acquisition” set Acquisition mode to “ran from PulseSpel”
- b. In the “experiment” chose “2P ESE set up”, phase cycle “ 2 step”. Make sure the set p0 and p1, attenuation and gain from previous steps.
- c. In pulseSPEL variable, type d0 and set it to 0.
- d. Type d1 and set it to 180 ns.
- e. Type srt, set it to 1000*1.02 ns
- f. Click the “run” in the left bottom corner and record the echo signal. Measure the time position of the point at the beginning of the echo signal which is located at approximately ¼ of the hight.
- g. After that, change the d0 to adjust the beginning of echo to this point.
- h. Measure the echo signal width at approximately ¼ height.
- i. Type pg in pulseSPEL variable and set this number for integration of echo signal, for example 74 for the nitroxide echo.
- j. Type d30 and set it to 4 ns.
- k. Change the dim2 to correct number for the echo decay. It should allow for seeing the whole echo decay to nearly zero.
- l. Make sure P0 and P1 are correct.
- m. Compile variable definitions and validate the program.
- n. Change the phase cycle to “2 steps”
- o. Click the “run” in the “left bottom corner”. You should get an ESR echo decay curve.

IMSERC Pulse Q-band Manual (V1.0)

2. Run "EDLOR optimization"
 - a. In the "experiment" chose "ELDOR pulse optimization", phase cycle "2 step".
 - b. Type d30 and set it to 2 ns.
 - c. Set P3 = 0 ns. Set af0 = 0 and af1 = 0, as0=0, as1=0 or 5. (0 is square pulse, 5 is chirp shape pulse).
 - d. Change the dim3 to a correct number for the nutation curve.
 - e. Compile variable definitions and validate the program.
 - f. Change the phase cycle to "2 steps"
 - g. Click the "run" in the "left bottom corner". You should get nutation curve and find the length of the pi pulse from it
3. Run "ELDOR" Experiment
 - a. In the "experiment" chose "4 pulse ELDOR Set up", phase cycle "2 step". Make sure you set p0, p1, p3, af1, attenuation and gain from previous steps.
 - b. Type d1 and set it to 180 ns. d2 = 1.2 us, d3=60 ns. In general, d2 should be chosen based on the echo decay curve: the signal should not completely decay and have some measurable intensity at double d2 time.
 - c. Type d30 and set it to 2 ns.
 - d. Set P3 = pi pulse from step 14 nutation curve. Set af0 = 0 and af1 = 80000
 - e. Change the magnetic field by adding 80/2.8=28.5 (conversion of the af1 value in MHz to Gauss).
 - f. Compile variable definitions and validate the program.
 - g. Change the phase cycle to "2 steps"
 - h. Click "run" in the "left bottom corner". You should get a refocused ELDOR echo followed some spurious signals from imperfect phase cycling.
 - i. Adjust d0 to the start of the signal. Typically d0 value will decrease by approximately 40 compared to step 13f.
 - j. Go to "Parameters" → "Patterns", click "Start". Go to SpecJet and click "Run".
 - k. Adjust the signal phase using SpecJet: maximize the green trace and minimize the yellow one.
 - l. Press "Stop" at patterns, press "Stop" on SpecJet. Change the phase cycle to "16 steps"
 - m. Click "run" in the "left bottom corner". You should get a single clean signal of the refocused ELDOR echo.
 - n. As previously, measure the width of the echo at ¼ height. Type "pg" and input this number as.
 - Click "run" in the "left bottom corner".
 - o. After the "4p pulse ELDOR set up optimization", switch to pulse "4P EDLOR 1D Acquisition"

IMSERC Pulse Q-band Manual (V1.0)

- p. Set d30 value (typically 4 or 8ns) according to the time length of your DEER trail and the desired number of points in it
- q. Make sure all the parameters in above steps are right, then calculate the sx number (calculated by the equation shown in the PulseSpel program) according to your d1, d2, d3 and d30 values.
- r. Change dim5 dimension with correct sx number.
- s. Compile variable definitions and validate the program.
- t. Change the phase cycle to "16 steps"
- u. Set the number of scans "n" according to your need.
- v. Click "run" in the "left bottom corner". You should get an ELDOR spectrum.

IMSERC Pulse Q-band Manual (V1.0)